
## **QUICK START GUIDE**

Review full instruction book prior to first use



## QUICK START GUIDE

waterpik waterflosser

Review full instruction book prior to first use

2 START WATER FLOSSING

TO AVOID MESS: PLACE THE TIP IN YOUR MOUTH BEFORE TURNING UNIT ON

Turn power and water on. Start on low pressure and increase to your liking. Close lips enough to prevent splashing, while still allowing water to flow from mouth into the sink.

4 TECHNIQUE

Aim water at the gumline at a 90° angle. Follow the gumline and pause briefly between teeth.

5

Sole Distributor (Hong Kong & Macau) Nudento International Limited 19/F., Yat Chau Building 262 Des Voeux Central Hong Kong Tel: 3176 0561

www.waterpik.com.hk

Model WP-861E